

Product Features:

- Clear fluid pathway facilitate priming during setup
- ☐ Ergonomic trigger flush device
- Modular transducer mounting plate system
- Transducer chip is integrated into the patient disposables

This device is indicated for clinical pressure monitoring use as transmission tubings.

Principle	Specification	Remarks
Utilizations of monocrystalline silicon's piezoresistive effect and electrical input/output conversion make it achieve precise measurement and monitoring of patient's blood pressure.	PC06-1-FRSC	
	PC06-2	
	PC06-3	



Purecath Medical (Shanghai) Co., Ltd. Add: 2A 6 Building, No. 328 Jinglian Road, 201108 Shanghai PEOPLE'S REPUBLIC OF CHINA

Tel.: +86-21-64348830 Fax: +86-21-64348829 Website:www.purecath.com